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Applications:

- For high-level instruments disinfectant
- Suitable for manual final disinfection of medical surgical and dental instruments as such as endoscopes
- Ready to use and doesn't need dilution
- Applicable for 50 times disinfection or 14 days of use

Material compatibility:

Suitable for Stainless steel, common plastics such as polycarbonates, PMMA and polystyrene, elastomers including polyurethane, silicone, Teflon, rubber, latex

Usage instructions:

- Infected instruments must be cleaned thoroughly before being disinfected with Denova.
- Immerse the instruments in a soaking tanks contains Denova (with the lid closed)
 - For strong disinfection: 5 minutes
 - For complete sterilization: 30 minutes
- Rinse instruments by clean water.

Physical Characteristics:


- Colorless Liquid
- Without any suspended particles
- Weak smell of acetic acid (vinager)
- Non-flammable

Corrosion Test:

Electrochemical Corrosion analysis and polarization plot

Electrode Type	Corrosion rate (mm/year)
Stainless steel 304	4.788×10^{-3}
Stainless steel 304	3.251×10^{-3}

* Razi Metallurgical Research Center - Tracking Number:12814

* Denova polarization plot

Storage:

- Store at room Temperature (4-30 °C)
- Keep out of reach of children
- Keep away from direct sunlight or high temperature
- Do not pour back used Denova into original container (clean Solution)
- Prevent skin and eye contact with Denova and wash hands or liquid contacted skin thoroughly with water after handling.

DeNova

Safe and effective peracetic acid based disinfectant for critical, semi-critical and non-critical medical, dental and surgical instruments.

Ingredients:

- Peracetic acid
- Acetic acid
- Hydrogen peroxide
- Adjuvants

Shelf Life:

- 24 months

Packaging:

- 5L container

Disinfection Efficacy:

Strains	Results After 5min contact (CFU/mL)
Sporicidal Bacillus Subtilis (ATCC6633)	-
Fungicidal Candida albicans (ATCC10231)	<5
Bactericidal Staphylococcus aureus (ATCC538) - Gram Positive	-
Pseudomonas aeruginosa (ATCC9027) - Gram Negative	-

* Tested at Pasteur institute of Iran- Control Number: 95-08-EC161